



Headquarters / Engineering research laboratory :
23 Gunpo Advanced Industry 1-ro(Bugok-dong), Gunpo-si, Gyeonggi-do
Tel +82-31-490-0800 Fax +82-31-490-0801

Yeongnam business office / Plant :
55 Gonghangap-gil, 85beon-gil, Gangseogu, Busan Metropolitan City
Tel +51-973-8518 Fax +51-973-8519

E-mail : info@gastron.com

www.gastron.com

GTD-2000TxW

Instruction Manual



Read in detail for correct use.

Gas & Flame Detection System

GASTRON



In case of a problem after purchasing the product,
please contact the address below.

· Address : 23 Gunpo Advanced Industry 1-ro, Gunpo-si, Gyeonggi-do
· Tel : 031-490-0800
· Fax : 031-490-0801
· URL : www.gastron.com
· e-mail : info@gastron.com

Thank you very much for purchasing a product from Gastron Co. Ltd.

Our Gastron Co.Ltd. is a company specialized in Gas Detector & Gas Monitoring System and have been recognized by many customers for the best quality and use convenience. We always seek to help our customers to find the product they need and we continuously research to develop gas detectors that satisfies our customers. From now on, you can solve all problems related to gas detectors with the products of Gastron. We Gastron will be responsible for your satisfaction.

This user manual describes operation and simple maintenance methods etc. for the gas detector. Please read carefully and store it in a safe place. Using it as a reference will help a lot when you have a question during use.

- For accurate operation of the gas detector, please perform an inspection and a calibration at least once per 6 months.
(* In reference to KOSHA GUIDE: P-135/6-2018 / 7.2 In-house inspection, section 2)
- For accurate operation of the gas detector, inspection and calibration using a reference gas before measurement is recommended.
- Failure to calibrate may result in malfunction of the equipment due to aging of the sensor.
- Only an individual specialized in gas detector may disassemble the device.
- Wire specification for the power cable must be decided in reference to the "Installation cable length" section.
- Please contact our company's technical support, e-mail, or website for inquiries related to inspection and calibration of the gas detector.

The product and manual are subject to change without a notice for the product's functional improvement and ease of use.

* KOSHA GUIDE : P-135/6-2018

Calibration must be performed at a frequency requested by the manufacturer and shall be performed quarterly when the calibration period is not specified.


| | | | |
|--|----|---|----|
| 1. Overview | 6 | 7.2. CALIBRATION MODE | 21 |
| 2. Configuration | 6 | 7.2.1. Zero Calibration | 21 |
| 3. Specification | 7 | 7.2.2. Span Calibration | 22 |
| 3.1. Basic Specifications | 7 | 7.3. ALARM Mode | 23 |
| 3.2. Mechanical Specifications | 7 | 8. Troubleshooting | 24 |
| 3.3. Electrical Specifications (Standard Type) | 8 | 8.1. Fault List | 24 |
| 3.4. Environmental Specifications | 8 | 8.2. Recovery List | 24 |
| 4. Name and description of each part | 9 | 9. Drawings and Dimensions | 26 |
| 4.1. Components | 9 | 9.1. Standard Type | 26 |
| 4.1. LAMP Function Description | 10 | 9.2. Warning light GTL-100 | 27 |
| 5. Installation | 11 | 9.3. When Connecting Warning Light | 28 |
| 5.1. Detachment of Housing Cover | 11 | 10. Precautions before Installation | 29 |
| 5.2. Main PCB Configuration | 12 | 10.1. Selecting a Place for Installation (Occupation Safety and Health Act Data) | 29 |
| 5.3. Terminal Configuration | 13 | 10.2. Selecting a Site for Installation (High-Pressure Gas Safety Control Act Data) | 29 |
| 5.3.1. Wiring for 4~20mA Source Operation Type | 14 | 10.3. Precautions during Installation | 29 |
| 5.3.2. Wiring for 4~20mA Sink Operation Type | 14 | 11. Revision History | 31 |
| 5.3.3. Wiring for 4~20mA 3Wire Sink Operation Type | 15 | | |
| 5.4. Installation Cable Length | 15 | | |
| 6. Detector Operation Flow | 17 | | |
| 6.1. Initial Operation Status (Power On) | 17 | | |
| 6.2. Measuring Mode | 17 | | |
| 6.3. Operation Flow | 18 | | |
| 6.4. Menu Configuration Table | 19 | | |
| 7. System Mode | 20 | | |
| 7.1. PROGRAM MODE | 20 | | |

GTD-2000TxW toxic gas detector has been developed to detect gas leaked from industrial sites and various toxic gases generated from factories, gas storages, and manufacturing processes that produce or use toxic gases and to prevent accidents in advance.

GTD-2000TxW toxic gas detector is installed in areas with toxic gas leak hazards and continuously monitors gas leak. It displays measurements on built-in displays (LCD or OLED) of the detector, delivers alarm status with visible and audible signals using a warning light, and provides DC 4~20 mA standard output signal. Also, for DC 4~20 mA standard output, output signal transmission length between detector and receiver can be connected up to 2,500 m. (When cable CVVS or CVVSB 1.5sq and higher is used.) GTD-2000TxW flammable gas detector must be used at height below 1,000 m above sea level.

2. Configuration

Body of GTD-2000TxW is made of Aluminum alloy and the gas sensor module is made of stainless steel. It consists of a complete explosion-proof enclosure (Ex d IIC T6). This product can be installed in areas with flammable gas leak and explosion hazards. It has built-in LCD on the detector to display gas leak status at installed site. It consists of display part that indicates measurements, terminal part that output measurements (DC 4~20mA) externally, and a PCB board. External configuration consists of sensor that detects gas leak, cable inlet, and high-pressure explosion proof small warning light made of SUS materials. It uses magnet-bar outside the main body of detector enabling calibration from the outside of the detector, thus, maintenance is convenient.



[Figure 1. GTD-2000TxW Overview]

3.1. Basic Specifications

| ITEMS | SPECIFICATION | |
|--------------------------|--|-------|
| Measuring Type | Diffusion | |
| Measuring Value Display | - LCD or OLED Display - Green, Red, Yellow LED - Sound Output(95dB) | |
| Measuring Method | - Electro-Chemical Cell - Heated-semiconductor Cell) - Photoionization detector(PID) | |
| Detectable Gas | Toxic Gas (Note 1) | |
| Measuring Range | Capable to display 000.0 ~ 9999 (Note 1) | |
| Accuracy | ≤ ±3% / Full Range | |
| Zero Drift | ≤2% / Full Range | |
| Response Time | Depends on Sensor Module. Refer to Sensor Specification or Contact in case for Special Gas. | |
| Approvals Classification | KCs: Ex d IIC T6,T5,T4, IP65 ATEX/IECEX: II 2 GEx d IIC Gb T4~T6 SIL2, MED, ABS, DNV | |
| Basic Interface | Analog 4-20mA current interface | |
| HART Interface | HART REV7 (Option) | |
| Option | HART Board | |
| Warranty | Transmitter & Lamp | 2Year |
| | Sensor | 1Year |

※ Note1. Refer to the measured gas list for measured gases and their ranges. Contact us for special gas.

3.2. Mechanical Specifications

| ITEMS | SPECIFICATION | |
|-------------------------|--|--------------------------|
| Explosion-Proof Type | Explosion-Proof Enclosure | |
| Dimension | - Transmitter: 136(W) × 166(H) × 110(D) mm - Lamp: 83.8(W) × 85(H) × 48(D) mm | |
| Weight including Sensor | App. 2.2kg | |
| Mounting Type | Wall mount | |
| Mounting Hole | Ø 7 ±0.1 | |
| Cable inlet | 3/4" PF (1/2" or 3/4" NPT) | |
| Body Material | Transmitter | Aluminum Alloy |
| | Sensor | Stainless Steel (STS316) |
| | Transmitter | PC |
| | Sensor | Stainless Steel (SUS304) |


3.3. Electrical Specifications (Standard Type)

| ITEMS | SPECIFICATION | |
|--|--|---------------------------|
| Input Voltage(Standard) ※ Customer supplied PSU must meet requirements IEC1010-1 and CE Marking requirements. | Absolute min: | 18V |
| | Nominal: | 24V |
| | Absolute max: | 31V |
| | Ripple maximum allowed: | 1V pk-pk |
| Wattage | Max. wattage: | 3.6W @+24 VDC |
| | Max. current: | 150mA @+24 VDC |
| Analog output Current | 0-20 mA (500 ohms max load) All readings ± 0.2mA Measured-value signal: 4mA (Zero) to 20mA (Full Scale) | |
| | Fault: | 0mA |
| | 0-100% LEL: | 4mA - 20mA |
| | 100-109%LEL: | 20mA - 21.4mA |
| | Over 110% LEL: | 22mA |
| | Maintenance: | 3mA |
| Analog output current ripple & noise max | ±20uA | |
| Wiring requirement | Power | CVVS or CVVSB with shield |
| | Analog | CVVS or CVVSB with shield |
| Signal Transmission Distance (Cable Connection Length) | Analog | 2500m |
| EMC Protection: | Complies with EN50270 | |

3.4. Environmental Specifications

| ITEMS | SPECIFICATION | |
|-----------------------|---------------|-------------------------------|
| Operation Temperature | Transmitter | -40 to 80 °C |
| | Sensor | Refer to Sensor Specification |
| Storage Temperature | Transmitter | -40 to 80 °C |
| | Sensor | Refer to Sensor Specification |
| Operation Humidity | Transmitter | 5 to 99% RH (Non-condensing) |
| | Sensor | Refer to Sensor Specification |
| Pressure Range | 90 to 110KPa | |
| Max. air velocity | 6m/s | |

4.1. Components



[Figure 2. GTD-2000TxW Components]

| No | ITEMS | SPECIFICATION |
|----|-----------------------|--|
| 1 | Housing Body | Protects PCB Board built in Sensor and Housing from external environmental change and shock. |
| 2 | Housing Cover | It is assembled with Detector Housing Body. Top surface is built with circular glass to enable monitoring of measurement displayed on LCD. |
| 3 | MAIN PCB | Amplifies fine outputs generated from Sensor Element to transmit a converted output in 4~20 mA DC standard. It sends data to display part. |
| 4 | Display PCB | Displays data sent from MAIN/ Terminal (Transmitter) PCB on LCD and displays power status with a power lamp. |
| 5 | Power/Signal Terminal | CN8 consists DC18V~31 power supply and DC 4~20 mA standard output connection terminal (VISO, +24V, mA, GND, ETH). |
| 6 | Sensor Terminal | CN4 is Sensor Connection Terminal. |
| 7 | Power LED | Lamp that turns on upon power supply. |
| 8 | Function Key | It is used to enter function setting mode from gas measuring mode by contacting magnet-bar for 2 sec or longer. It is also used to save changed data from function setting mode. |

| No | ITEMS | SPECIFICATION |
|----|-----------------------|--|
| 9 | Reset Key | To cancel or return to the previous status during Parameter setting, use the Magnet-Bar and touch once. (Each touch returns to the previous status by one unit.) |
| 10 | FUNC key | During conversion of mode or number, use the Magnet-Bar and touch once. Each touch converts or increases displayed value by one unit. |
| 11 | ↓(DOWN)Key | During conversion of mode or number, use the Magnet-Bar and touch once. Each touch converts or increases displayed value by one unit. |
| 12 | External Earth Ground | - It must be grounded to outside of detector for protection from external noise or strong electric field. - Use a conductor that is 4 mm or longer when coupling ground line. |
| 13 | Mount Hole | Hole to fix the gas detector on external wall or other installation sites. |
| 14 | Sensor | It is a site that detects actual gas leak. It converts the amount of gas leak into electrical signal and transmits to the AMPPCB. |
| 15 | Conduit Connection | It is supplied for inlet of power supply and measurement output signal for the detector during installation. For cable inlet, PF or NPT 1/2", 3/4" are prepared. |
| 16 | Internal Ground | - It must be grounded to inside of detector for protection from external noise or strong electric field. - Use a conductor that is 4 mm or longer when coupling ground line. |
| 17 | FRONT CASE | It is an explosion-proof case made of explosion-proof plastic and emits LED light externally through transparent window. |
| 18 | CASE RING | It serves as a cover to prevent deviation of front case. |
| 19 | SEXANGLE BOLT | It serves to prevent deviation after case ring assembly and is fixed by mood bolt. |
| 20 | CASE BODY | It is the main body that protects internal buzzer and other parts by body. |
| 21 | BUZZER | It receives signal upon an event of error and generates warning sound. |
| 22 | FIXED RING | It prevents deviation of sinter filter upon external shock. |
| 23 | FLAME ARRESTER | It blocks flow in of dust, debris, water, gas and delivers warning sound from the buzzer to outside. |
| 24 | ARREST CAP | It protects sinter filter from external shock. |
| 25 | LED COVER | It is a joint for LED PCB ASS'Y and covered with moulding. |
| 26 | LED PCB ASS'Y | It is a multi-color LED that displays status. |
| 27 | NUT | It serves to fix parts to the compliment part. |
| 28 | SHAFT BODY | It serves to fix parts to the compliment part. |
| 29 | CABLE WIRE | It supplies power and transmits data. |

[Table 1. GTD-2000TxW 구성요소 설명]


4.2. LAMP 기능설명

- In normal status, green LED lights on and buzzer does not ring.
- When Alarm-1 signal is received, red LED flickers in 1 sec interval and buzzer rings in 1 sec interval.
- When Alarm-2 signal is received, red LED lights on and buzzer rings continuously.
- When Trouble signal is received, yellow LED lights on and buzzer does not ring.

It is prohibited for an individual, other than an approved user or a technician responsible for installation and repair from the head office, to install a gas leak detector on site or open the cover of the installed gas leak detector and manipulate it. This may cause serious loss of life and property from fire, explosion, and etc. In addition, please check whether there is any remaining explosive gas or combustible material in the surroundings. Power must be turned off before performing work.


5.1. Detachment of Housing Cover

- <Warning - Do not open when electrical current is flowing>
- Turn the slotted set screw (M4 x 1ea) fixing the cover part of main body 3~4 turns counter clockwise (ccw) using a hex wrench (M2) then turn the cover of gas leak detector ccw to detach the cover.
When the cover is detached, LCD appears.



[Figure 3. Slotted Set Screw]


- After detaching the cover, disassemble the display parts as below.
 - ① Push in left and right fixing rings located on front side of LCD at the same time.
 - ② While pushing, pull the display parts towards the front to detach from gas detector body.
 - ③ After detaching the display parts, the Main PCB is installed at the bottom part of the detector body.



[Figure 4. Display Part Detaching Method]

5.2. Main PCB Configuration

■ After detaching the display parts, the Main PCB terminal layout appears as shown in the figure below.



[Figure 5. Main PCB Terminal Layout]


| No | NAME | DESCRIPTION |
|----|-------------|--|
| 1 | CN9 | Power & Output Signal Terminal |
| 2 | J1 | 4~20mA Source / Sink Selection Jumper (ON: Source Type, OFF: Sink Type) |
| 3 | CN1 | Display LCD Connector |
| 4 | CN8 | Sensor Connector |
| 5 | CN6 | Program download Connector |
| 6 | D1 | Status LED (Flashes in 1 sec interval during normal operation) |
| 7 | OC1,OC2,OC3 | HART Option Board Connector |
| 8 | CN4 | Warning light (GTL-100) Interface Connector |

[Table 2. Main PCB Key Part Description]

■ Using OC1, OC2, and OC3 shown in the layout above, HART option board can be attached and 3ØScrew holes located at top left of HART option board are used for fixing.

5.3. Terminal Configuration

- <Warning - Turn off power before connecting power terminal>
- After disassembling display parts, there is a terminal block in the Main PCB as shown in the figure below. Holding it with hands and pulling towards ceiling detaches it from the Main PCB.
- Loosen 5 terminal fixing screws located at top part of detached terminal block CN9 (VIS, +, mA, -, ETH) Connector by turning counter-clockwise using a Ø driver. Connect DC 18~24 V power to +, and - then connect signal cable to mA. Tighten 5 terminal fixing screws clockwise so that terminal does not leave the track then insert Main PCB as the same condition before disassembly.



[Figure 6. CN9 Terminal Configuration]


| NO | PCB SILK | PIN NAME | DESCRIPTION | |
|----|----------|----------|------------------------------------|-----------------------------------|
| | | | 4~20mA Source Drive (J1 Jumper ON) | 4~20mA Sink Drive (J1 Jumper OFF) |
| 1 | VISO | VIS | N.C | 4~20mA Sink In(+) |
| 2 | +24V | + | +24V / POWER (+) | |
| 3 | mA | mA | 4~20mA Source Out | 4~20mA Sink Out(-) |
| 4 | GND | - | GND / POWER (-) | |
| 5 | ETH | ET | EARTH | |

[Table 3. CN9 Terminal Detailed Description]

- Use CVVS or CVVSB 2.0sq⊥ Shield Cable for terminal configuration.
- To connect 4 Pin terminal from the existing old GTD-2000TxW model, connect terminals in reference to pin #2, which is +24V.

5.3.1. Wiring for 4~20mA Source Operation Type


- Connect 4-20 mA signal terminal at PLC side to 'mA' of GTD-2000Tx. GND terminal is used in common with power. Then, turn on the J1 jumper.
- HART Communicator can only be used in models using HART Option board.



[Figure 7. 4-20mA Source Configuration]

5.3.2. Wiring for 4~20mA Sink Operation Type


- Connect (+) and (-) terminals for 4-20 mA sink output at PLC side to VISO terminal and 'mA' terminal, respectively. Then, turn off the J1 jumper.
- HART Communicator can only be used in models using HART Option board.



[Figure 8. 4-20mA Sink Configuration]

5.3.3. Wiring for 4~20mA 3Wire Sink Operation Type


- Connect (+) and (-) terminals for 4-20 mA sink output at PLC side to VISO terminal and power (24V DC) (-) terminal, respectively. Connect 'mA' terminal of GTD-2000TxW to 'GND' terminal. Then, turn off the J1 jumper.



[Figure 9. 4-20mA 3Wire Sink Configuration]

5.4. Installation Cable Length

- The maximum length between GTD-2000TxW and power supply is decided by wire specification.
- Max. Installation Length = $V_{MAXDROP} \div I_{MAX} \div WIRER/m \div 2$
 - $V_{MAXDROP}$: Maximum Power Loop Voltage Drop (=Power Supply voltage - min operating voltage)
 - I_{MAX} : Max. Current of GTD-2000TxW
 - $WIRER/m$: The resistance of the wire (ohms/meter value available in wire manufacturer's specification data sheet)
- Example of installation lengths using 24 V power supply and 16 AWG is as follows.
 - GTD-2000TxW minimum operating voltage = 18 Vdc
 - $V_{MAXDROP} = 24 - 18 = 6V$
 - $I_{MAX} = 0.15A (150mA)$
 - $6 \div 0.15 \div 0.01318 \div 2 = 1,517.451m \approx 1,517m$



[Figure 10. Calculation of GTD-2000TxW Installation Cable Length]

■ Power cable installation for each cable type is as shown in the table below.

| AWG | mm ² | COPPER RESISTANCE(ohms/m) | METERS |
|-----|-----------------|---------------------------|--------|
| 12 | 3.31 | 0.00521 | 3838 |
| 14 | 2.08 | 0.00828 | 2415 |
| 16 | 1.31 | 0.01318 | 1517 |
| 18 | 0.82 | 0.02095 | 954 |
| 20 | 0.518 | 0.0333 | 600 |

[Table 4. GTD- 2000TxW Power Cable Installation Length]

6.1. Initial Operation Status (Power On)

■ After wiring to power terminal at the top of Main PCB board then supply power, the following contents are displayed on LCD. Approx. 30 min of stabilization of time is needed from the initial supply of operation power and it starts to operate normally after sufficient stabilization.

| | |
|--------------------------|---|
| GTD-2000 [V X.XX] | - When power turns on, model name and product firmware version are displayed on LCD (OLED). |
| SELF TEST [>>>0179] | - Self test is performed for 3 min (30 min for O2 set-up), '>' symbol in the second row shows the progress. |


6.2. Measuring Mode

■ After power on, when there is no error from "SELF TEST", it automatically enters Measuring Mode.

| | |
|------------------------|---|
| [%LEL] [0.0] | - In normal operation mode, it runs as below. - In the first row, measured gas name and unit are displayed in 1 sec interval. In the second row, the current measurement is displayed. |
| [COMB.] [0.0] | - In the current screen, touching the magnet-bar on 'Func' key for 2 sec or longer changes to setting mode. ※ During Modbus network, '*' is displayed on the left of the first row. ※ ENG. When the mode is on, '<' is displayed on the left of the second row. |
| [%LEL] [30] | - When alarm1 or 2 occurs, it runs as below. (It only runs only when alarm used setting in alarm mode is on.) - The first row runs the same as in normal status. In the second row, alarm message and gas measurement are displayed in 1 sec interval. |
| [%LEL] [ALARM1] | - When GTL100 explosion-proof warning light is installed, upon an event of alarm1, red LED and buzzer flicker and run in 1 sec interval. Upon an event of alarm2, they continuously run without flickering. - If alarm latch is on, the alarm continuously run until it is released using the reset key. |
| [%LEL] [OVER] | - When gas measurement larger than 10% of set high scale value, "OVER" displays in 1 sec interval. At this time, 4~20 mA runs as 22 mA. |
| [%LEL] [UNDER] | - When gas measurement smaller than -10%, "UNDER" displays in 1 sec interval and 4~20 mA runs below 2 mA. ※ The following function only runs when UNDER is on. |
| [FAULT8] SEN HUNT | - When trouble occurs in the device, it displays fault number and message. - At this time, 4~20 mA runs below 2 mA. - Fault #8 on the left is a mode that is displayed when there is repeated hunting error from the sensor. |

6.3. Operation Flow

- After power on, it passes self-diagnostic process then enters Measuring Mode. Here, by operating front keys, you can go to internal System Mode.
- Timeout for Level1 and Level2 are 10 sec. It is set to 1 h for Level2 Calibration and Test Mode.
- When "RESET" key contacts at Program Mode Screen, it returns to Measuring Mode. When "RESET" key contacts at each Program Setting Screen, it returns to the parent step.



[Figure 11. GTD-2000TxW Mode Configuration]

- Operation keys for detector system mode are defined as below.

| ITEM | NAME | DESCRIPTION |
|-------|--------------|---|
| FUNC | Function key | Detector Mode Setting Entry Function (Insert Magnetic-bar for 2 sec or longer in measuring mode) Level2 Next stage entering function and setting saving function |
| RESET | Reset key | Move to the previous stage from the level entered. |
| ↑ | Up key | Next level mode that is configured in LEVEL1 and Change in Level2 setting Plus |
| ↓ | Down key | Next level mode that is configured in LEVEL1 and Change in Level2 setting Minus |

6.4. Menu Configuration Table

| LEVEL1 | LEVEL2 | | DEFAULT |
|-------------------------------------|---------------------------------------|-----------------------------|------------|
| | NAME | PARAMETER | |
| PROGRAM MODE | GAS TYPE (Gas Type) | [DEFIN./USER] | DEFIN. |
| | GAS SEL (Gas Select) | Built-in gas name selection | COMB. |
| | UNIT SEL (Unit Select) | %/%LEL/PPM/PPB | %LEL |
| | D-POINT (Decimal Point) | 0.100/1.00/10.0/100 | 100 |
| | HIGH SCL (High Scale) | 1~10000 | 100 |
| | PASSWORD | 00~99 | 00 |
| CALIBRA. MODE (Calibration Mode) | CALIBRA. [ZERO] (Calibration Zero) | ZERO CAL [NO,YES] | [NO] |
| | | ZERO GAS [0] | |
| | | ZERO >>>> [SUCCESS / FAIL] | |
| | | CAL. DATA [0 / FAIL] | |
| | CALIBRA. [SPAN] (Calibration Span) | SPAN CAL [NO, YES] | [NO] |
| | | SPAN GAS [0] | |
| | | SPAN SET [50 / FAIL] | 50 |
| | | SPAN >>>> [SUCCESS / FAIL] | |
| | | CAL. DATA [0] | |
| | ALARM MODE | ALM USED (Alarm Used) | [OFF / ON] |
| AL LATCH (Alarm Latch) | | [OFF / ON] | OFF |
| A1 LEVEL (Alarm-1 Level) | | [1 ~ Full SCL] | 20 |
| A1 TYPE (Alarm-1 Type) | | [INC / DEC] | INC |
| A1 DBAND (Alarm-1 Dead Band) | | [0.0 ~ 10.0%F.S] | 0 |
| A1 TIME (Alarm-1 Time) | | [0~ 60]SEC | 1 |
| A2 LEVEL (Alarm-2 Level) | | [1 ~ Full SCL] | 40 |
| A2 TYPE (Alarm-2 Type) | | [INC / DEC] | INC |
| A2 DBAND (Alarm-2 Dead Band) | | [0.0 ~ 10.0%F.S] | 0 |
| A2 TIME (Alarm-2 Time) | | [0~ 60]SEC | 1 |
| SENSOR MODE | Factory Setting | | |
| MAINTEN.MODE | Factory Setting | | |
| DEVICE MODE | Factory Setting | | |
| VERSION MODE | Factory Setting | | |
| TEST MODE | Factory Setting | | |

7.1. PROGRAM MODE

| | |
|-------------------------------|---|
| PASSWORD [**] | - Contacting "FUNC" key with the Magnet-bar for 2 sec or longer in Measuring Mode enters Password mode. - After setting Password using "↑" or "↓" key, contact "FUNC" key. |
| PROGRAM MODE | - If password is correct, it enters Program mode. - By contacting "↑" or "↓" key, mode changes in defined order. (PROGRAM -> CALIABRA. -> ALARM -> SENSOR -> MAINTEN. -> DEVICE -> VERSION -> TEST) |
| GAS TYPE [DEFIN,] | - Mode that sets gas name type and contacting "↑" or "↓" key changes gas type. (DEFIN. / USER) - DEFIN. Is acronym for Define and is selected when using built-in setting of gas name. - USER is selected when the user directly sets a gas name. - When a desired gas type is displayed, contact "FUNC" key to set the gas name and enter Program Mode. |
| GAS SEL [COMB.] | - Mode that sets gas name and contacting "↑" or "↓" key changes gas name. - When gas type is set to DEFIN., built-in gas names can be selected for use. When it is set to USER, the user shall set gas name using 5 character for use. Characters allowed for use are number, capital Alphabet, space, and period. When position is moved by setting characters, a cursor in black box shape appears as shown in the figure on left. |
| USER GAS [SER] | |
| UNIT SEL [%LEL] | - Mode that sets gas measuring unit and contacting "↑" or "↓" key changes gas measuring unit. (% / %LEL / PPM / PPB). - When a desired gas measuring unit is displayed, contact "FUNC" key to set the unit and enter Program Mode. |
| D-POINT [100] | - Mode that changes decimal place and contacting "↑" or "↓" key changes decimal place. (0.100/1.00/10.0/100) - When a desired decimal place is displayed, contact "FUNC" key to set the decimal place and enter the next Program Mode. |
| HIGH SCL [100] | - Mode that sets high scale value that should be displayed for full range. Contacting "↑" or "↓" key increases or decreases the scale value. (1 ~ 10000) - When a desired high scale is displayed, contact "FUNC" key to set the high scale value and enter the next Program Mode. |
| PASSWORD [00] | - It is password setting mode. This password is checked when entering Program Mode. - Set using "↑" or "↓" key then contact "FUNC" key to set the password and enter the next Program Mode. |

7.2. CALIBRATION MODE

- Due to characteristics of the gas detector, minimum 30min of stabilization time is required and maintenance condition may change depending on site condition.

7.2.1. Zero Calibration

| | |
|-----------------------------------|--|
| PASSWORD [**] | - Contacting "FUNC" key with the Magnet-bar for 2 sec or longer in Measuring Mode enters Password mode. - After setting Password using "↑" or "↓" key, contact "FUNC" key. |
| CALIBRA. MODE | - Contact "↑" or "↓" key to select "Calibration Mode". - Contact "FUNC" key when "CALIBRA. MODE" is displayed to enter Calibration Mode. - Contact "RESET" key to return to Measuring Mode. |
| CALIBRA. [ZERO] | - When "CALIBRATION MODE" is displayed, contacting "FUNC" key selects Zero Calibration. - Contact "↑" or "↓" key to achieve [ZERO] then contact "FUNC" key to enter Zero Calibration. |
| ZERO CAL [NO] | - Contact "↑" or "↓" key to achieve [YES] then contact "FUNC" key to perform Zero Calibration. |
| ZERO GAS [0] | - Using a calibration tool, inject clean air or 100% nitrogen into the sensor at a flow rate of 500 mL/min for 1 min. Contact "FUNC" key when measurement is stabilized to automatically perform Zero Calibration. |
| ZERO >> [SUCESS] | - When zero calibration is successful, "ZERO SUCCESS" is displayed for 2 sec on LCD display then it changes to "CALIBRATION DATA" Mode. - When zero calibration is not successful, "ZERO FAIL" is displayed for 2 sec and it changes to "Calibration Data Mode". |
| ZERO >> [FAIL] | - ZERO FAIL occurs when sensor input value exceeds 70% of total ADC input range. |
| CAL.DATA [0] | - It is a mode to display measurement after calibration and checks whether calibration is successful after performing auto calibration. - When it failed, it displays "FAIL" and the current measurement in 1 sec interval. - Contact "RESET" key to return to "Calibration mode". |

7.2.2. Span Calibration

| | |
|------------------------------------|--|
| CALIBRA. MODE | <ul style="list-style-type: none"> - Contact "↑" or "↓" key to select "Calibration Mode". - Contact "FUNC" key when "CALIBRA. MODE" is displayed to enter Calibration Mode. - Contact "RESET" key to return to Measuring Mode. |
| CALIBRA. [SPAN] | <ul style="list-style-type: none"> - Contact "↑" or "↓" key to achieve [SPAN] then contact "FUNC" key to enter Span Calibration mode. |
| SPAN CAL [NO] | <ul style="list-style-type: none"> - Contact "↑" or "↓" key to achieve [YES] then contact "FUNC" key to perform Span Calibration. |
| SPAN GAS [50] | <ul style="list-style-type: none"> - Using a calibration tool, inject the standard gas to the sensor at a flow rate of 500 mL/min for 90 sec. After the measurement is stabilized, press "FUNC" key to enter the next mode. - When Hold function in Maintenance Mode is on, the current maximum SPAN gas value is held and displayed in the first row. |
| [50] [50] | <ul style="list-style-type: none"> - The current measurement is displayed in the second row. |
| SPAN SET [50] | <ul style="list-style-type: none"> - It is a mode to set standard gas value. When there is no "Fail" message, contact "↑" or "↓" key to set a value. - When the inject standard gas value is abnormal, it displays Fail message. Fail message is as follows. Fail message and span set value are displayed alternatively. |
| SPAN SET LOW FAIL | <ul style="list-style-type: none"> ① "LOW FAIL" occurs when the difference between injected gas value and zero calibration value is below 1%. |
| SPAN SET HIGH FAIL | <ul style="list-style-type: none"> ② "HIGH FAIL" occurs when injected gas value is above 95% of total ADC input range. - "RING FAIL" occurs when injected gas value runs in the current set span value, it exceeds 95% of ADC input range. The above FAIL messages are automatically released if the status is normal after adjusting SPAN value. |
| SPAN >> [SUCCESS] | <ul style="list-style-type: none"> - After completing standard gas value setting, contacting "FUNC" key automatically runs Span Calibration. When it succeeds, "SPAN SUCCESS" is displayed on LCD display for 2 sec then it changes to "CAL DATA" Mode. |
| SPAN >> [FAIL] | <ul style="list-style-type: none"> - When span calibration is not successful, "SPAN FAIL" is displayed for 2 sec and it changes to "Cal. Data Mode". |
| CAL.DATA [0] | <ul style="list-style-type: none"> - It is a mode to display measurement after calibration and checks whether calibration is successful after performing auto calibration. - Contact "RESET" key to return to "Calibration mode". |

7.3. ALARM MODE

| | |
|----------------------------|--|
| PASSWORD [**] | <ul style="list-style-type: none"> - Contacting "FUNC" key with the Magnet-bar for 2 s or longer in Measuring Mode enters Password mode. - After setting Password using "↑" or "↓" key, contact "FUNC" key. |
| ALARM MODE | <ul style="list-style-type: none"> - Contact "↑" or "↓" key to select "Alarm Mode". - Contact "FUNC" key when "ALARM MODE" is displayed to enter Alarm setting mode. - Contact "RESET" key to return to Measuring Mode. |
| ALM USED [ON] | <ul style="list-style-type: none"> - It is a setting where Alarm Mode Setting is turned on or off. - Contact "↑" or "↓" key to change ON/OFF status. When it is ON, Alarm function and Alarm Mode setting can be performed. - GTL-100 explosion-proof warning light can only be used when it is ON. - Contacting "FUNC" key enters the next mode. |
| AL LATCH [ON] | <ul style="list-style-type: none"> - It is a mode that sets Reset method after running of Alarm1. Contact "↑" or "↓" key to change "ON" and "OFF". - "OFF" sets to automatically reset the alarm. "ON" set to reset the alarm only when Reset key is pressed. - Contact "FUNC" key when a desired mode is displayed to set then enter the next alarm setting mode. |
| A1 LEVEL [20] | <ul style="list-style-type: none"> - It is a mode that sets Alarm1 level. Contact "↑" or "↓" key to increase or decrease Alarm1 level, respectively. - When a desired alarm 1 threshold is displayed, press "FUNC" KEY to set it and enter the next item. |
| A1 TYPE [INC] | <ul style="list-style-type: none"> - It is a mode that sets operational direction of Alarm2. Contact "↑" or "↓" key to display "INC" or "DEC". - "INC" mode operates when the value is equal or larger than set alarm threshold. "DEC" mode operates when the value is equal or less than set alarm threshold. |
| A1 TYPE [DEC] | <ul style="list-style-type: none"> - Contact "FUNC" key when a desired mode is displayed to set then enter the next mode. |
| A1 DBAND [0]%F | <ul style="list-style-type: none"> - It is a mode that sets Dead band value for Alarm1 operation. Use "↑" or "↓" key to set a value. - It is a function where Alarm1 operates at values above the sum of Alarm1 level and dead band values and releases at values below Dead band value subtracted from Alarm1 level. - When a desired alarm 1 dead band is displayed, press "FUNC" KEY to set it and enter the next item. |
| A1 TIME 0 SEC | <ul style="list-style-type: none"> - It is a function to prevent instantaneous malfunction of detector due to external shock and noise other than from normal operation and time can be set in a range between 0~60 sec. - For Alarm1 dwell time setting, press "↑" or "↓" key to increase or decrease in unit of 1 sec, respectively. - When a desired Alarm 1 dwell time is displayed, press "FUNC" KEY to set it and enter the next item. - Ex.) Alarm threshold value: 20% LEL / Delay time: When it is at 5 sec, Alarm triggers when the measured value is above the set value based on 20%LEL for 5 sec or longer. When it goes down below the set value within 5 sec, alarm is not triggered. |
| A2 LEVEL [40] | <ul style="list-style-type: none"> - It is a mode that sets Alarm2 level. Contact "↑" or "↓" key to increase or decrease Alarm2 level, respectively. - When a desired alarm 2 threshold is displayed, press "FUNC" KEY to set it and enter the next item. |

| | |
|---------------------------|--|
| A2 TYPE [INC] | - It is a mode that sets operational direction of Alarm2. Contact "↑" or "↓" key to display "INC" or "DEC". - "INC" mode operates when the value is equal or larger than set alarm threshold. "DEC" mode operates when the value is equal or less than set alarm threshold. |
| A2 TYPE [DEC] | - Contact "FUNC" key when a desired mode is displayed to set then enter the next mode. |
| A2 DBAND [0]%F | - It is a mode that sets Dead band value for Alarm2 operation. Use "↑" or "↓" key to set a value. - It is a function where Alarm2 operates at values above the sum of Alarm2 level and dead band values and releases at values below Dead band value subtracted from Alarm2 level. - When a desired alarm 2 dead band is displayed, press "FUNC" KEY to set it and enter the next item. |
| A2 TIME 0 SEC | - It is a function to prevent instantaneous malfunction of detector due to external shock and noise other than from normal operation and time can be set in a range between 0~60 sec. - For Alarm2 dwell time setting, press "↑" or "↓" key to increase or decrease in unit of 1 sec, respectively. - When a desired Alarm 1 dwell time is displayed, press "FUNC" KEY to set it and enter the next item. - Ex.) Alarm threshold value: 40% LEL / Delay time: When it is at 5 sec, Alarm triggers when the measured value is above the set value based on 40%LEL for 5 sec or longer. When it goes down below the set value within 5 sec, alarm is not triggered. |

8.1. Fault List

| FAULT MESSAGE | DESCRIPTION & CONDITION | CAUSE |
|----------------------|---|---|
| FAULT2 "SEN HIGH" | Sensor output is above ADC max. value. | Defective sensor module or transmitter board ADC |
| FAULT3 "SEN LOW" | Toxic sensor output is below ADC min. value. | Defective sensor module or transmitter board ADC |
| FAULT4 "EROM ERR" | Defective Transmitter EEPROM Checksum | Defective Transmitter board EEPROM |
| FAULT5 "+24V LOW" | It occurs when 24 V main input power is inputted with voltage below 17V. | Input power below 17 V or defective Transmitter ADC |
| FAULT6 "ADC FAIL" | Defective Transmitter ADC | Defective Transmitter ADC |
| FAULT7 "H/W REV" | H/W Version Error | Defective MPU inside Transmitter |
| FAULT8 "SEN HUNT" | When there is repeated hunting error in gas measurement, sensor output status may be poor or check for input power noise. | |


[Table 5. Fault List]

8.2. Recovery List

| No | CAUSE | SOLUTION |
|----|------------------------------------|---|
| 1 | Defective Sensor Module | Change Sensor |
| 2 | Defective Transmitter Board ADC | Change Transmitter Board |
| 3 | Defective Transmitter Board EEPROM | 1) Perform Factory Initialization then correct parameter and re-calibrate 2) Change sensor unit when the same problem occurs again |
| 4 | Defective MPU inside Transmitter | Change Transmitter Main Board |
| 5 | Poor Sensor Output Status | Change Sensor |
| 6 | Check Input Power Noise | Check External Input Voltage Noise Status |


[Table 6. Recovery List]

9.1. Standard Type




[Figure 12. GTD-2000TxW Standard Type Drawing]

9.2. Warning Light GTL-100



[Figure 13. GTL-100 Drawing]

9.3. When Connecting Warning Light



[Figure 14. GTD-2000TxW Warning Light Connection Drawing]

10.1. Selecting a Place for Installation (Occupational Health & Safety Act Data)

A gas leak detector alarm shall be installed in the following places.

- Around chemical equipment and accessories that have concerns of gas leak. This includes compressors, valves, reactors, pipe joints, etc. installed inside and outside of a building that handle combustible and toxic materials.
- Places that are easier for gases to stay such as areas around manufacturing facilities with ignition sources like heating furnace, etc,
- Areas around equipment for filling combustible and toxic materials.
- Substations, panel rooms, control rooms, and etc. located within explosive area.
- Other areas that are easier for gases to stay.

10.2. Selecting a Site for Installation (High-Pressure Gas Safety Control Act Data)

Gas detector of gas leak detector alarm must be installed as close to the areas with concerns of gas leakage as possible. However, for areas where direct gas leakage is not expected but are easier for leaked gas to stay, the detector must be installed at the point 1 of the following.

- Gas leak detector alarm installed outside a building shall be installed at points where gas is likely to stay in consideration to wind direction, wind speed, specific gravity of gas, etc.
- Gas leak detector alarm installed inside a building shall be installed near the floor when the specific gravity of gas is heavier than air and near ventilation of ceiling when it is lighter than air.
- Alarm for gas leak detector alarm must be installed at sites where the gas detector is installed and workers are present.


10.3. Precautions during Installation

Avoid areas with electrical barriers such as rain water, etc. It is recommended to be installed in areas that are easier to work in since regular maintenance is needed. Avoid areas with vibration or shock since they can affect output values. Sensor part must be installed towards the direction of gravity.


- This equipment has explosion-proof construction for internal pressure and belongs to GROUP II for gas and vapor in general work sites and chemical plants. It can be used in ZONE 1 (ONE) and ZONE 2 (TWO) hazardous sites.
- Allowable temperature is 85 C or below, which corresponds to T6.
- Surrounding temperature shall be in a range of -40 to 80 °C (for Main body. For sensor, refer to each sensor manual.)
- Installation Height: 1,000 M below sea level
- Relative Humidity: 5% ~ 99% (Non-condensing)
- Installation Site: Indoor and Outdoor
- Explosion Ignition Group for Target Gas or Vapor: Ex d IIC T6

- During wiring work, use explosion-proof cable gland at cable inlet or tightly seal cable conduit during metal cable wiring construction to prevent spread of flames in case of explosion or movement of gas, etc. through the cable conduit within 50 mm. All materials including materials used for sealing of unused inlets must have safety certificates!
- When connecting the equipment with cable, screw thread must be tightened 5 threads or more.
- Work in conditions satisfying other [Standards for Selection, Installation, and Maintenance, etc. of Explosion-proof Electric Machine and Equipment Wiring, etc. at Work Site]

| VERSION | CONTENTS | DATE |
|---------|---------------------------|------------|
| 0.1 | * Manual Initial Revision | 2017.09.29 |



[Figure 15. High-Pressure Packing Type]



[Figure 16. Y Sealing Compound]